


Handicare 4000 Curved Stairlift

SEAT MEASUREMENTS







Please note, the diagrams above provide a legend for measurements and do not represent the exact details of any specific seat.









		Simplicity	Style	Smart*	Perch
A	Swivel clearance	29" (735 mm)	29" (735 mm)	27.5 - 29.5" (700 - 750 mm)	-
	with short armrest	27.5" (700 mm)	27.5" (700 mm)	_	_
B	Seat depth	15.25" (390 mm)	15.25" (390 mm)	14.25" - 16.75" (360 - 425 mm)	14.25" - 16.75" (360 - 425 mm)
C	Rail depth from wall	5" (124 mm)	5" (124 mm)	5" (124 mm)	5" (124 mm)
D	Fold-up depth	16" (405 mm)	16" (405 mm)	16" (405 mm)	16" (405 mm)
(3)	Fold-out depth	27.4" (695 mm)	27.4" (695 mm)	27.4" (695 mm)	27.4" (695 mm)
	with small footrest	25" (636 mm)	25" (636 mm)	25" (636 mm)	25" (636 mm)
	with offset footrest	25.75" (655 mm)	25.75" (655 mm)	25.75" (655 mm)	25.75" (655 mm)
•	Width between armrests	17" (430 mm)	17" (430 mm)	16.25 - 20" (410 - 510 mm)	16.25 - 20" (410 - 510 mm)
G	Backrest height	15.75" (400 mm)	15.75" (400 mm)	18" (455 mm)	18" (455 mm)
•	Footrest to seat	18.5 - 21.75" (470 - 550 mm)	18.5 - 21.75" (470 - 550 mm)	19.25 - 23.25" (490 - 590 mm)	24.5 - 29.5" (625 - 750 mm)
	with higher footrest	16.5 - 19.75" (420 - 500 mm)	16.5 - 19.75" (420 - 500 mm)	17.25 - 21.25" (440 - 540 mm)	_
0	Footrest height	2.75" (70 mm)	2.75" (70 mm)	2.75" (70 mm)	2.75" (70 mm)

^{*}The Smart seat is also available in UV-protected fabric for outdoor use (in slate color).



To find a dealer near you, call us at 1.888.637.8155 handicareusa.com | handicare.ca

Handicare is a division of Savaria Corporation